



## Tobra Specimen Basket

Sterile, safe, continuous, and efficient collection of specimen.

### KEY FEATURES



Unique mesh basket allows for confident collection of specimen

Mesh basket allows for continuous filtration of specimen that does not interrupt surgery



Safe removable basket for easy transport of specimen

Offering a **user friendly** one step design for easy set up and efficient performance. Designed **small and compact** for in-line collection.



Specimen may be transported in the device cup by connecting supplied tubing to both barbs

Additional tubing is supplied to facilitate Specimen Basket connection with devices such as aspirators and shavers to suction canisters.



Interested in easy to use bone collection?  
Tobra Part #BB-100 (5 units per box)

CALL US TODAY 866.777.1505